

Key

Microscope Worksheet:

Calculating Magnification. Converting Measurements, Estimating cell size, Calculating Field of View & Magnification Problems

1. Calculate total magnification: Ocular x Objective

| Ocular | Objective | Total Magnification |
|--------|-----------|---------------------|
| 10X | 4X | 40X |
| 15X | 10X | 150X |
| 5X | 12X | 60X |
| 10X | 10X | 100X |
| 10X | 40X | 400X |

2. What are the possible magnifications of a microscope with an ocular marked 10X and objectives marked 5X, 15X, 30X and 60X?
50X, 150X, 300X, 600X

3. Convert the following measurements:

- a. $9.2 \text{ mm} = \underline{9200} \text{ } \mu\text{m}$
- b. $5900 \text{ } \mu\text{m} = \underline{5.9} \text{ mm}$
- c. $0.083 \text{ mm} = \underline{83} \text{ } \mu\text{m}$
- d. $61000 \text{ } \mu\text{m} = \underline{61} \text{ mm}$

4. Estimating cell size: Divide the field of view by the number of cells that occupy the diameter.

- a. The field of view is 2500um. If a cell takes up 1/5 of the field of view, how long is the cell?

$$\text{Estimated size} = 2500\text{um}/5 \text{ cells} = 500\text{um}$$

- b. A student counts 50 cells across the diameter of the field of view, and there are 70 rows of cells. If the diameter of the field of view is 3500 μm , what is the length and width of the cells?

$$\text{Length} = 3500\text{um}/50 = 70\text{um}$$

5. Calculate the field of view (Field Diameter):

- a. Low power: $4X = \frac{4500 \text{ } \mu\text{m}}{1000} = \underline{4.5} \text{ mm}$
- b. Medium power: $10X = \frac{1800 \text{ } \mu\text{m}}{1000} = \underline{1.8} \text{ mm}$
- c. High Power: $40X = \frac{450 \text{ } \mu\text{m}}{1000} = \underline{0.45} \text{ mm}$

8. The diagrams below represent what can be seen through an imaginary microscope. You have been given information along with the diagrams. Using the information provided, calculate the size of the objects viewed.

a) This triangle is being viewed under medium power.


low power objective = 4X

medium power objective = 10X

high power objective = 40X

eyepiece = 5X

low power FOV = 4.2 mm = 4200 μm



a) Calculate Size of Medium F. of V.

$$\frac{\text{L.P. Mag}}{\text{M.P. Mag}} \times \frac{\text{L.P. Field of View}}{1} = \frac{4}{10} \times \frac{4200}{1} = \frac{16800}{10} = 1680 \mu\text{m}$$

b) Calculate Size of Object (Longest way)

$$\text{Size} = \frac{\text{FofV}}{\# \text{ across}}$$

$$\text{If you think 3 fit} \rightarrow \frac{1680 \mu\text{m}}{3} = 560 \mu\text{m}$$

$$\text{If you think 4 fit} \rightarrow \frac{1680 \mu\text{m}}{4}$$

$$\leftarrow \text{Anything Between} \rightarrow = 420 \mu\text{m}$$

b) This diamond is being viewed under high power.


low power objective = 5X

medium power objective = 40X

high power objective = 100X

eyepiece = 10X

low power FOV = 3.5 mm



$$\text{Size} = 58.3 \mu\text{m} \leftrightarrow 43.8 \mu\text{m}$$